

価格表：

製品番号	製品名	説明	テスト数	希望小売価格
1-100-10010	Q-Fusion Leukemia Screening Kit	30 fusion genes	10	¥277,000
1-100-10020	Q-Fusion Leukemia Screening Kit	30 fusion genes	20	¥553,000
1-100-20010	AML Q-Fusion Screening Kit	16 fusion genes with 72 breakpoints	10	¥236,000
1-100-20020	AML Q-Fusion Screening Kit	16 fusion genes with 72 breakpoints	20	¥472,000
1-100-30010	ALL Q-Fusion Screening Kit	16 fusion genes with 71 breakpoints	10	¥236,000
1-100-30020	ALL Q-Fusion Screening Kit	16 fusion genes with 71 breakpoints	20	¥472,000

機器適合表：

	Detection Channels	BioRad CFX96 FAM/HEX / ROX/CY5	ABI7500 FAM/HEX/ ROX/CY5	Agilent Mx3005p FAM/HEX/ ROX/CY5	Roche LC480 FAM/HEX/CY5	Qiagen Rotorgene FAM/HEX/ ROX/CY5	ABI StepOne FAM/HEX/ ROX
Leukemia Q30	FAM/HEX/ ROX/CY5	+	+	+	X	X	X
AML/ALL	FAM/HEX/ ROX	+	+	+	X	X	+
BCR-ABL1 One Tube Quant	FAM/HEX/ CY5	+	+	+	+	+	X
Lung Cancer ALK/ROS1/RET Testing Panel	FAM/HEX/ CY5	+	+	+	+	+	X
EML4-ALK	HEX/CY5	+	+	+	+	+	+
Lung Cancer EGFR/KRAS/BRAF Mutation Panel	FAM/HEX/ CY5	+	+	+	+	+	+
Q-MRSA	FAM/HEX/ ROX/CY5	+	+	+	X	X	X
HPV Genotyping	FAM/HEX/ ROX/CY5	+	+	+	X	X	X
TB	FAM/TET	+	+	+	+	+	+

価格はいずれも税別です。
本製品は、研究用試薬としてのみご利用いただけます。

お問合せ： info2@rikengenesys.jp
営業・企画部：03-5759-6042